

Qiagen Gel Extraction - quick protocol

Prepared by Ms Alex Aitken

1. Run a TAE agarose gel as normal – 1% maximum is best.
2. Excise the gel slice with clean blades and place in eppendorff tube.
3. Add QG buffer to the top of the tube – add at least 1ml; if it does not look like 1ml will fit, cut the gel slice into 2 pieces and put into 2 tubes and process both tubes as individual samples.
4. Incubate the gel slice in QG buffer at 50°C for approx 10 minutes, vortex every couple of minutes.
5. Ensure that the gel has completely dissolved.
6. Pipette 750µl into one of the columns.
7. Centrifuge the column for a few secs (typically up to top speed and then stop is sufficient).
8. Discard the flow through and add the remaining sample to the column.
9. Repeat #7.
10. Discard the flow through.
11. Add 750µl of PE buffer (with alcohol added as per instructions on bottle) to the column and leave to stand at RT for 1 minute.
12. Repeat #7.
13. Discard the flow through.
14. Centrifuge the column for 1 minute at top speed (this ensures that all the alcohol is removed from the column).
15. Place the column in a fresh tube.
16. Add 25-50µl TE (or use dH₂O directly onto the disc at the base of the column (do not touch the disc)).
17. Stand at RT for 5 minutes.
18. Centrifuge at top speed for 1 minute.
19. Discard column.
20. Store eluted DNA at -20 °C.